

Multi-parameter Patient Monitor



PM-MP-1-1

Features & Benefits:

- * 12.1" color TFT display, high resolution: 800×600dpi;
- * Separated para board: ECG, NIBP, SpO2 board;
- * Full touch-screen enables intuitive operation by clicking one the specific parameters or the waveform in real time, standard configuration includes the navigation knob.
- * Real time S-T segment analysis, and ARR analysis;
- * SpO2 Pitch Tone Variation and Drug Dose Calculation;
- * Parameter waveform and characters color selectable
- * Multi-display: Standard, Large font, Trend Coexis, OxyCRG dynamic
- * Standby mode for Energy saving and suspend monitoring;
- * 7-lead ECG waveforms display in phase;
- * Capture dynamic waveforms;
- * Built-in rechargeable lithium battery, and battery volume checking,
- * AC power checking
- * Efficient resistance to interference of defibrillator and HF knife;
- * SINNOR F-6 SpO2 technology, VS NELLCOR accuracy;
- * Networking capacity and TCP/IP networking platform enable abundant future upgrade
- * Time delay switch-off function keep the monitor data restart monitor;

Multi-parameter Patient Monitor



- * Complete and flexible mounting solution for transportation and in-patient monitoring
- * Standard parameters: 5-lead ECG, SpO2(F-5),NIBP, TEMP,RESP,PR,HR.
- Optional: Touch Screen, Thermal Recorder, 3/12-lead ECG, 2/4-IBP,4-TEMP,2-SpO2, Nellcor/Masimo SpO2,BIS, Phaisein/ Respironics Multi-gas, ICG/CO, FHR, FM, TOCO.

| Performance Specifications | |
|--|--|
| Display | 12.1"color TFT ;Rolling and refreshing waveform display |
| Resolution | 800×600 |
| Multi display selectable including | Large-font display/Trend coexist display OxyCRG dynamic view display./Bed-to-bed view display |
| Trace | 9 waveforms(7 ECG, 1 SPO2 and 1 RESP) |
| Sweep speed | 12.5mm/s,25mm/s,50mm/s |
| Indicator | Power/battery indicator light QRS beep and alarm sound |
| Battery | Rechargeable lead acid cell, 12v/4AH Max.24 hour for charging, 4 hours for continue working |
| Trend/Parameter graphic and tabular trends | 5s/piece,8 hours; 1min/piece,168 hours(24hour×7days) 5min/piece,1000 hours. |
| Storage | NIBP: 1000 groups Alarm: 200 groups Full disclosure waveforms: 3600S SD card external storage |
| Alarm | User-adjustable High, Medium and Low limits 3-level Audible and visual alarm |
| Recorder | Build-in, thermal array Plethysmogram caveform:2 channels Record mode: manual, on alarm, time-defined Recording width: 50mm Printing speed: 50mm/s |
| Recording type | Frozen waveform record NIBP recall record Trend table record Alarm record Fixed -time record |

Multi-parameter Patient monitor



| ECG | |
|----------------------|--|
| Lead mode: | 5 –lead(R,L,F,N,C) |
| Lead selection : | I,II,III,avR,avL,avF,V |
| Waveform: | 3 and 7 channel selectable |
| Gain selection: | 0.5mm/mv,1mm/mv,2mm/mv |
| Sweep speed: | 12.5mm/s;25mm/s;50mm/s |
| Heart rate range: | Adult: 15~300bpm; Neonate/pediatric:15~350bpm |
| Accuracy | +1bpm or +1%,whichever is greater |
| Resolution | 1bpm |
| Filter | surgery mode:1~20Hz monitor model:0.5~40Hz Diagnostic mode:0.05~130Hz |
| Scaling signal | 1mv,+3% |
| Protection | withstand 4000VAC/50 voltage isolation against electrosurgical interference and defibrillation |
| Alarm range | 15~350bpm |
| S-T segment detecton | |
| Measurement range | 2.0mV~+2.0mV |
| Alarm Range | -2.0mV~ +2.0mV |
| Accuracy | -0.8mV ~+0.8Mv |
| Error | ±0.02Mv |
| Arrhythmia analysis | YES |
| SPO2 | |
| Measurement range | 0~100% |
| Resolution | 1% |
| Accuracy | ±2%(70-100%);0-69% unspecified |
| Alarm range | 0~100% |
| range | 20~300bpm |
| Resolution | 1bpm |
| Error | ±1bpm or ±2%,whichever is greater |

Multi-parameter Patient monitor



| NIBP | |
|-----------------------------|---|
| Method | Digital Automatic oscillometric |
| Operation mode | Manual/Automatic/continuous |
| Auto measurement time | Adjustable(1~480min) |
| Measurement Unit | mmHg/Kpa selectable |
| Measurement types | Systolic, Diastolic, Mean |
| Measurement tange | |
| Range of Systolic pressure | Adult: 40~270mmHg Pediatric: 40~220mmHg Neonate: 40~135mmHg |
| Range of mean pressure | Adult: 20~235mmHg Pediatric: 20~165mmHg Neonate: 20~110mmHg |
| Range of diastolic pressure | Adult: 10~215mmHg Pediatric: 10~150mmHg Neonate :10~100mmHg |
| Over-pressure protection | Double safety protection |
| Resolution | 1 mmHg |
| Alarm | Systolic.Diastolic,Mean |
| Respiration | |
| Method | Thoracic impedance |
| | Adult:7~120rpm; Neonate/Pediatric:7~150rpm |
| Apnea alarm | YES,10~40s |
| Resolution | 1rpm |
| Accuracy | ±2rpm |
| Temperature | |
| Compatible probe | YSI or CYF |
| Measurement range | 5~50°C |
| Resolution | 0.1°C |
| Accuracy | ±0.1°C |
| Refreshing time | about 1 |
| Average measuring time | <10s |

Multi-parameter Patient monitor



| FHR | |
|---|--|
| Transducer | Multi-crystal, Pulsed Doppler |
| Measurement range | 50~210 BPM |
| Working frequency | 1 MHz |
| Strength | <5mW/cm[] |
| Signal processing | special DSP system and modern recognition. |
| Resolution | 1BPM |
| Accuracy | ±1BPM |
| Alarm Range | High: 160,170,180,190 BPM, Low: 100,110,120 BPM |
| FM | |
| Manual button marking, automatic FM identifying function | |
| TOCO Measuring | |
| Transducer | External pressure transducer |
| Measurement range | 0~100 units |
| Resolution | 1rpm |
| Accuracy | ±2 rpm |
| IBP | |
| Channel | 2 channels |
| Range | -50-300mmHg |
| Resolution | 1mmHg |
| Accuracy | ±4mmHg(±4%) |
| Unit | mmHg, Kpa |
| Converter Sensitivity | 5V/V/mmHg |
| Transducer sites | ART/PA/CVP/LAP/RAP/ICP |
| EtCO2(Sidertream CO2) | |
| Measurement range | 0~99mmHg |
| Accuracy | ±2mmHg (0~40mmHg) |
| Sampling Range | 100ml/min |
| Sampling rate accuracy | 15% |
| Respiration rate | 0~120rpm |
| Respiration accuracy | ±2rpm (0~70rpm);+5rpm (>70rpm) |
| Respiration time | <240msec (10% to 90%) |
| Delay time | <2s |

Multi-parameter Patient monitor



| EtCO2 (mainstream CO2) | |
|---------------------------------|--|
| Method | Infrared Spectrum |
| Range | 0.0-10% (0~76%) |
| Resolution | 1mmHg (0.1%) |
| Accuracy | <5%(±4.0 mmHg);Or<10%(of Readings) |
| Recorder | |
| Build-in, thermal array | |
| Plethysmogram caveform | 2 channels |
| Record mode | manual, on alarm, time-defined |
| Recording width | 50mm |
| Printing speed | 50mm/s |
| Recording type | Frozen waveform record NIBP recall record Trend table record Alarm record Fixed -time record |
| Satety | |
| Safety level | Class I, type CF |
| Dimention and Weight | |
| Dimention | 28*42*48cm |
| G.Weight | <6.5KS |
| Operation Environment | |
| Temperature | Working 0~+40°C Transportation and storage -20~+60°C |
| Humidity | working≤85% Transportation and storage≤93% |
| Power | AC 100-240,50/60Hz |
| Patient range | Neonate,pediatric,and adult patients |
| Accessories supplies | |
| (1) 5 lead ECG cable | |
| (2) 1 spo2 probes | |
| (3) 1 NIBP prbe | |
| (4) 1 temp probe | |
| (5)1 ground connecting liner | |
| (6) Chest electrode (10pcs/set) | |